Technical Sheet December 2017

Wrongwoods

Richard Woods & Sebastian Wrong



Type: Storage

1090


2600



355

1500

510



Design: 2007

Designers: Richard Woods & Sebastian Wrong

Dimensions

Night Table: $L585mm \times D355mm \times H595mm$ Chest of Drawers: $L1090mm \times D510mm \times H735mm$ Low Cabinet: L1500mm × D405mm × H645mm

L1000mm × D340mm × H1460mm Room Divider:

Long Low Cabinet: L2600mm × D405mm × H645mm

Weight

16kg 70kg 60kg 140kg

70kg

Materials

Plywood, timber, paint, clear lacquer

Notes: Artisan material and process, variations occur

Colour Options

Night Table:

Pink with blue, white with black, white with green, yellow with red

Chest of Drawers:

Pink with blue, white with black, white with green, yellow with red

Room Divider:

Light grey with dark grey, white with yellow

Long Low Cabinet:

Pink with blue, white with black, white with green, yellow with red

Low Cabinet:

Pink with blue, white with black, white with green, yellow with red, light grey with dark grey

L1070mm × D380mm × H1500mm

Colour References

Pink with Blue NCS S1040R, NCS S1050 R90B, NCS S0560 Y20R

Yellow with Red NCS 1080 Y80R, NCS 2070 Y60R, NCS 0580Y20R, NCS2070 Y20R Green NCS S2070 G10Y, NCS1075 G40Y, NCS S2060 G70 Y, NCS S2070 G60Y

Grey NCS S77502-B, NCS S5502-B, NCS S2002-B, NCS S0500-N Black with White RAL 9005 - RAL 9010

Packaging Dimensions		Weight	Packaging Type	Notes: part assembled, 2 people
Night Table:	L630mm × D400mm × H500mm	18kg	Cardboard	required for handling and assembling
Chest of Drawers:	L1180mm × D600mm × H800mm	105kg	Crate	
Low Cabinet:	L1550 mm × D450mm × H550mm	91ka	Crate	

Crate

Crate

182kg

115kg

Care & Maintenance

Room Divider:

This piece is not suitable for outdoor use.

Do not expose to direct sunlight or high temperatures.

Long Low Cabinet: L2650mm × D450mm × H550mm

Avoid contact with heat sources.

Please use a slightly damp microfibre cloth, avoiding any abrasives or solvents.

Please remember to remove any jewellery when cleaning. In order to guarantee the stability of the article, periodically check the legs are tightened.

